



Epley manoeuvre for treating left-sided benign paroxysmal positional vertigo (BPPV) [NB1]


Clinician instructions




1. Put a pillow a quarter of the way down the examination couch. Ask patient to sit in the middle of the couch.




3. Ask patient to look up and tip the head back slightly, then turn the head 45 degrees to the left.




5. Wait in this position for 1 minute.




7. Turn head a further 90 degrees to right, while telling patient to roll onto right side. Wait in this position for 1 minute.




2. Swing patient's legs onto the couch so they are straight out in front.



 Keep head in this position and lie patient down with head over top edge of pillow—pillow is at shoulder level and head is about 30 degrees below the horizontal.



6. Turn head 90 degrees to right. Wait in this position for 1 minute.



 Sit patient up slowly. Tell patient not to lie down flat again until bedtime. Warn them that they may have mild nausea or dizziness or feel unsteady for several hours after a successful treatment.

NB1: Steps 1 to 4 are the same as a Hallpike manoeuvre.